

System Description - T-101

Document No.	SYS-T101-001	Revision	01
System	Solution Preparation Tank T-101	Status	Approved - Training
Project	Otrivin Nasal Spray Facility - Synthetic Case	Purpose	RAG/C&Q Training

Training disclaimer: This is a fictional, simplified document set created solely for education. It does not describe or reproduce the actual manufacturing process, specifications, or controlled records for Otrivin.

1. Purpose

T-101 is the solution preparation vessel in this synthetic xylometazoline nasal-spray training facility. The vessel receives purified water and formulated ingredients, provides controlled mixing, and transfers prepared solution to the downstream filtration/holding step.

2. Major equipment

- 1,000 L nominal working-volume vessel T-101
- Variable-speed agitator AG-101
- Temperature transmitter TT-101
- Level transmitter LT-101
- CIP spray device
- Product outlet valve XV-104
- PLC/HMI recipe control

3. C&Q; boundary

The exercise boundary starts at the T-101 inlet connections and ends at the product outlet valve XV-104. Utilities and downstream filling equipment are outside this package except where interfaces must be verified.