

Vendor Manual Extract - T-101

Document No.	VM-T101-001	Revision	01
System	Solution Preparation Tank T-101	Status	Approved - Training
Project	Otrivin Nasal Spray Facility - Synthetic Case	Purpose	RAG/C&Q Training

Training disclaimer: This is a fictional, simplified document set created solely for education. It does not describe or reproduce the actual manufacturing process, specifications, or controlled records for Otrivin.

1. Operation

Before use, verify the vessel is released for operation, the outlet valve is in the required position, and the agitator is available. Charge liquid before starting the agitator.

2. Agitator

Permitted operating range is 50-300 rpm. Operation below minimum speed may not provide the intended mixing performance.

3. Cleaning

The vessel is equipped with a fixed spray device intended for automated CIP. Cleaning cycle parameters are controlled by the site automation system.

4. Maintenance

Inspect agitator seal and product-contact gaskets at defined preventive-maintenance intervals.